

HAKUZO

360 Degree Toothbrush



360 degree bristle



To remove plaque
with capillary action

For tongue cleaning

Product Features

With the use of 24,000 0.06 mm ultra thin bristles (20x of regular toothbrush), plaque can be removed efficiently with capillary action.

Capillary action is a phenomenon that liquid is sucked up between thin tubes by the action of surface tension of the liquid.



Instructions Do not rotate the brush. Brush the teeth like regular toothbrush.

Front teeth



Backside of the front teeth



Cuspid



The upper and lower front teeth



Back teeth



The upper and lower back teeth



Gum care




Tongue cleaning



Shaft Material	PP	Bristles Type	Soft
Bristles Material	Nylon	Heat Resistant Temperature	90°C

Specifications

Product Code	Specifications	Qty/box	Qty/carton	JAN Code
3162360	Brush dia. 15mm White	10 brushes	50boxes	 4 958286 718187



Design Your Voice
Hakuzo

Add: 2-4-9, Tokui-cho, Chuo-ku,
 Osaka 540-0025, Japan

TEL: 06-6942-0451

FAX: 06-6945-0539

www.hakuzo.co.jp/eng/index.php