

HAKUZO 360 Degree Toothbrush



贅沢毛量
20
TIMES
MORE
BRISTALS!

極やわらかブラシ

360°
Design



360 degree bristals


To remove plaque
with capillary action

For tongue cleaning

Product Features

With the use of 24,000 0.06 mm ultra thin bristles (20x of regular toothbrush), plaque can be removed efficiently with capillary action.

Shaft Material	PP	Bristles Type	Soft
Bristles Material	Nylon	Heat Resistant Temperature	90°C



- The 360-degree bristles allow for cleaning from any angle.
- High-density bristles can remove plaque effectively.



- Employing ultra-fine bristles with rounded tips to prevent gum damage.



- The compact brush head can reach deep into the mouth.
- It cleans uneven surfaces effortlessly.



Capillary action is a phenomenon that liquid is sucked up between thin tubes by the action of surface tension of the liquid.

Instructions Do not rotate the brush. Brush the teeth like regular toothbrush.

Front teeth



Backside of the front teeth



Cuspid



The upper and lower front teeth



Back teeth



The upper and lower back teeth



Gum care



Tongue cleaning



Specifications

Product Code	Size and Color	Qty/Box	Qty/Carton	JAN Code
3118067	Brush dia. 15mm, White	30 pcs	16 boxes	4 958286 720388

Made in Japan

Scan the QR code
to watch the
instructional video



Design Your Voice

Hakuzo

Add: 2-4-9, Tokui-cho, Chuo-ku,
Osaka 540-0025, Japan

TEL: 06-6942-0451

FAX: 06-6945-0539

www.hakuzo.co.jp/eng/index.php