

Toothbrush specifically designed for assisted oral care

# ComfortFit Brush



We are here to support you with your oral care concerns!

- ✓ Oral care for individuals with trismus
- ✓ Oral care for individuals with sensitive teeth

The brush is designed to closely conform to the surface of teeth and gums

Slightly angled handle conforms to the teeth and gums closely.



Easily slips into the oral cavity, even for individuals with difficulty opening their mouths.

Soft


Ultra Soft



# Easy Access to Oral Cavity


※The illustration is for reference only.

## Regular Toothbrush



Regular toothbrushes can cause discomfort during oral care by pulling at the cheeks.

## ComfortFit Brush



Slightly angled handle allows easy access to the oral cavity and conforms to the teeth and gums closely.

## Ultra soft and soft bristles

2 types of bristles are available, designed to meet diverse needs and preferences.

For individuals with sensitive teeth or those experiencing gum inflammation.




### Ultra Soft



Bristle Color: White / Handle Colors: Orange and yellow

### Soft



Bristle Color: Blue / Handle Colors: Blue and green

Can remove plaque effectively, suitable for individuals with healthy gums and teeth.



## Small Brush Head

The small brush head ensures easy access to hard-to-reach areas in oral cavity, allowing for thorough cleaning of each tooth.

## Hexagonal Handle

The hexagonal handle design is optimized for a comfortable pen-like grip, enhancing both control and cleaning efficiency

Handle Material	PP
Bristle Material	Nylon
Heat resistance temperature	80°C



Product code	Specifications	Qty/box	Qty/carton	JAN Code
3118054	Ultra Soft, 2 pcs/pack	12 packs	12 boxes	4 958286 720890
3118055	Soft, 2 pcs/pack	12 packs	12 boxes	4 958286 720906



Design Your Voice  
**Hakuzo**

Add: 2-4-9, Tokui-cho, Chuo-ku,  
Osaka 540-0025, Japan

TEL: 06-6942-0451

FAX: 06-6945-0539

[www.hakuzo.co.jp/eng/index.php](http://www.hakuzo.co.jp/eng/index.php)