

《50 pieces pack X 2 alcohol swabs》



Wiper Twin Pack Series

Wiper Twin Pack Series is a 50 pieces pack X 2 containers packaging.

Wiper Twin Pack Series is low in quantity per packaging and low in costs.

Elepha Cotton Eco E

Ethanol swab for disinfection



This alcohol swab is impregnated with ethanol approved by Japanese Pharmacopoeia.

【Contents and amount】

Medical absorbent cotton	9.3g
Japanese Pharmacopoeia approved ethanol (76.9~81.4vol%)	55mL

※Amount shown is for 1 container (50 pcs)

Elepha Cotton Eco ISO

70vol% isopropanol swab



This alcohol swab is impregnated with 70vol% isopropanol.

【Contents and amount】

Medical absorbent cotton	9.3g
Isopropanol (70vol%)	55mL
distilled water	q. s.

【Product Features】

【Container/Packaging】

2 containers can be separated easily just by tearing.

【Economical】

Small packaging is usually referred as high costs, however Twin Pack helps to lower down the costs with its relatively small packaging.

【Convenience】

50 pcs packaging is small in packaging and easy to carry.

【Hassle free】

Save time in preparing solutions and also swabs, directly saving costs.




【How to use】

- Separate the 2 containers, open the cover and tear open the aluminium foil.
- Close the cap properly after using.

Hakuzo Medical Corporation


How to use an exclusive Hakuzo Twin Pack fixed container

STEP 1




Open the container cap (Please do not tear open the aluminium foil).

STEP 2




Keep the container cap opened. Slide the container in slowly and make sure that it is completely slid in.

STEP 3




Close the container cap.

STEP 4




Peel off the seal out from the adhesive that is attached at the back of the exclusive Hakuzo Twin Pack fixed container.

STEP 5




The process is completed once the exclusive Hakuzo Twin Pack fixed container is fixed. (Only one hand is needed to take out the swabs).



When removing the exclusive Hakuzo Twin Pack fixed container. Pull it from the back when trying to remove the container.

The stickiness of the adhesive may reduce with the increase in the number of use, please wash it by water and the stickiness will recover.



HAKUZO

Wiper Twin Pack Series

《50 pieces pack X 2 alcohol swabs》

⚠ Attention

❌ Instructions (Please follow the instructions strictly)

Please avoid these areas

- Wound.
- Eye surrounding, mucous membrane (Eye, mouth, nose and etc).

👤 Consultation related

1. For the following candidates, please consult a doctor or a pharmacist.
 - A patient who is receiving medical treatment.
 - A patient who is allergy to medications, food and etc.
2. If the following symptoms occur, please stop using immediately and bring along this leaflet to consult a doctor, or pharmacist.

Related Part	Symptoms
Skin	Rashes, rubefaction, itching

Other Attention Effect on clinical test

The usage of this product may cause allergy before the actual allergy test injection takes place, directly affecting the results of the allergy test.

Other functions Can be used for disinfection on fingers, skin, and also medical devices.

< How to use >



- Wipe on the desired part.

< Instructions and dosage >

- For external use only
- Avoid contact with eyes. In case of contact with eyes occurs, wash eyes with water immediately. Please consult a doctor immediately if situation worsens.
- When applied on a wide area or used for a long period of time, inhaling may harm.
- When applied on the same area repeatedly, it may cause skin to get dried.
- Adult's supervision is needed when being used by a child.
- Coagulin may occur, please wash the medical equipment accordingly when coagulin occurs.

< Storage >

- Avoid direct sunlight and store in a cool place.
- Keep out of children.
- Do not refill with other packaging.
(might cause misuse and deterioration in quality)
- It may cause color fading on clothing when contacted.
- Please use up as soon as possible once the seal is opened.
- Please keep out of fire.

Product Name	Product Code	Specifications	Packaging	Qty per case	JAN code
Elepha Cotton E	2600101	4cm × 4cm	50pcs × 2containers/set	40 sets	 4 958286 330488
Elephone Cotton ISO	2600104	4cm × 4cm	50pcs × 2containers/set	40 sets	 4 958286 330518



Hakuzo Medical Corporation

4-9, 2-chome, Tokuicho, Chuo-ku, Osaka
540-0025, Japan
Tel: 81-6-6809-4319 Fax: 81-6-6945-0539
E-mail: hmc-trade@hakuzo.co.jp
<http://www.hakuzo.co.jp/eng/>