



HAKUZO

Powder Free

PLASTIC GLOVES EL

Hygienic powder free plastic gloves

100 PCS



Product Features


Hygienic powder free gloves. No risk of contamination caused by powder splashes.

Avoid skin rash caused by powder. No messy residue on hands or clothes.


Hakuzo Medical Corporation

INSTRUCTIONS

How to Put on



How to Take off



Material

- PVC

Cautions

- This product may cause skin allergies like rash or itchy skin. Stop using immediately when feeling discomfort.
- Keep hands away from high temperature heat sources when wearing this product.
- Avoid contact with concentrated solution and solvent like gasoline and benzene.
- Avoid contact with oil or food containing oil.
- Be careful not to damage the gloves when handling sharp tools.
- Avoid direct sunlight and hot and wet environment. Store in a clean place.

Product Code	Size	Qty/box	Qty/carton	JAN Code
3024101	S	100 pcs	20 boxes	4 958286 706269
3024102	M	100 pcs	20 boxes	4 958286 706276
3024103	L	100 pcs	20 boxes	4 958286 706283



Hakuzo PPE Holders

Customized holders that can be used with PPE products are sold separately.



Design Your Voice

Hakuzo

Add: 2-4-9, Tokui-Cho, Chuo-ku, Osaka 540-0025
 TEL: +816-6809-4319
 FAX: +816-6942-0496
 www.hakuzo.co.jp