

Hakuzo Fine Guard Face Shields

Disposable face shields for
protecting face from splashes.

The glasses frame fits nose bridge firmly
and can be worn together with eyeglasses
and a N95 mask.

Can be worn together
with eyeglasses.

Can protect users from
splashes from above.




with
Anti-Fog
Coating

Regular

with
Anti-Fog
&
Anti-Reflective
Coating

AR

Glasses frames are sold separately



10 each

Blue/Yellow/Pink



Hakuzo Fine Guard AR 10 PCS

Light and thin shield with anti-reflective and anti-fog coating that can be worn for a long period of time. Suitable for computer and EHR users.



PET 0.125mm (with protective film)

Hakuzo Fine Guard Regular 50 PCS

Environment-friendly recycled PET is being used. With anti-fog coating, the shield is suitable for short-term daily use.



PET 0.2mm



Assembly Method

Scan the QR code to watch the instructional videos.

Regular



AR



Product Code	Item	Qty/Box	Qty/Carton	JAN Code
3002799	Set Regular	Frame x 1 (Blue) Shields x 5	50 boxes	4 958286 717982
3002800	Set AR	Frame x 1 (Blue) Shields x 5	50 boxes	4 958286 717883
3002801	Frames (Blue/Yellow Pink)	30 pcs (10 each)	12 boxes	4 958286 717890
3002802	Shields Regular	50 pcs	15 boxes	4 958286 717906
3002803	Shields AR	10 pcs	100 boxes	4 958286 717913

Made in Japan



Design Your Voice

Hakuzo

Add: 2-4-9, Tokui-cho, Chuo-ku, Osaka 540-0025, Japan

TEL: 06-6942-0451

FAX: 06-6945-0539

www.hakuzo.co.jp/eng/index.php