

A disposable surgical mask  
with high barrier filter



Design Your Voice  
**Hakuzo**

# Hakuzo Omega Mask **Plus**



Pink  
**50** pcs

Disposable Surgical Mask



Minimising discomfort on ears with  
**soft elastic bands**

### Product Features

- Comes with metal nose strip and unique design on both sides help fits the face well.
- Elastic bands are made of latex-free materials.

### [MATERIALS]

Main body, filter: Polypropylene  
Elastic bands: Polyurethane  
metal nose strip part: Polyethylene

### [Test data]

Bacterial filtration efficiency ( average droplet size of 3.0 microns ) BFE 99%  
Particulate filtration efficiency ( average particle size of 0.1 micron) PFE 99%  
Fluid Resistance to synthetic blood : Above 21.3kPa  
Delta P (mm H2O/cm2): < 4.0

Test carried out by: KAKEN TEST CENTER General Incorporated Foundation

Product code	Type	Size	Packaging	Qty/case	JAN code
3901403	Pink	Regular 17.5cx x 9cm	50pcs/box	40 boxes	 4 958286 706245



Design Your Voice  
**Hakuzo**

Add: 2-4-9, Tokuicho, Chuo-ku  
Osaka 540-0025, Japan  
TEL: 0081-6-6942-0451  
FAX: 0081-6-6945-0539  
[www.hakuzo.co.jp/eng/index.php](http://www.hakuzo.co.jp/eng/index.php)

A disposable surgical mask  
with high barrier filter



Design Your Voice  
**Hakuzo**

# Hakuzo Omega Mask **Plus**

**50** pcs

Disposable Surgical Mask



Minimising discomfort on ears with  
**soft elastic bands**

## Product Features

- Comes with metal nose strip and unique design on both sides help fits the face well.
- Elastic bands are made of latex-free materials.

### [MATERIALS]

Main body, filter: Polypropylene  
Elastic bands: Polyurethane  
metal nose strip part: Polyethylene

### [Test data]

Bacterial filtration efficiency ( average droplet size of 3.0 microns ) BFE 99%  
Particulate filtration efficiency ( average particle size of 0.1 micron ) PFE 99%  
Fluid Resistance to synthetic blood : Above 21.3kPa  
Delta P (mm H2O/cm2): < 4.0

Test carried out by: KAKEN TEST CENTER General Incorporated Foundation

Product code	Type	Size	Packaging	Qty/case	JAN code
3901421	White	Regular 17.5cx x 9cm	50pcs/box	40 boxes	4 958286 711386
3901423	Blue	Regular 17.5cx x 9cm	50pcs/box	40 boxes	4 958286 711393



Design Your Voice  
**Hakuzo**

Add: 2-4-9, Tokuicho, Chuo-ku  
Osaka 540-0025, Japan  
TEL: 0081-6-6942-0451  
FAX: 0081-6-6945-0539  
[www.hakuzo.co.jp/eng/index.php](http://www.hakuzo.co.jp/eng/index.php)