

HAKUZO INFECTION CONTROL

W SURFACE DISINFECTANT WIPES BLOCK NEO

Large Sheet

EPA certified quaternary ammonium salt is being used.

Neutral Detergent



Press the button to open the lid.



60 count



30 count



250 count Refill Pack



250 count Bucket

PURPOSES

To disinfect and clean environmental surfaces such as table, door knob, bed rail, etc.



Table



Door knob



Bed rail

Size 150 × 300mm

Wipe material Rayon nonwoven fabric

Ingredients Quaternary ammonium salt, water

Hakuzo Medical Corporation

Product Feature No. 1

Solution approved by US EPA is being used. Effective in getting rid of a wide range of bacteria and germs.

Product Feature No. 2

Neutral and alcohol free solution barely cause effect on surfaces wiped. With low volatility, wide surface area can be wiped with 1 sheet.

Product Feature No. 3 60-count pack can be used with a customized holder. Convenient and user-friendly.

Screw fixable holder



Holder with hooks



Holder + pipe clamp (φ15mm~32mm)



Can be fixed on nursing trolley, bed rail, etc.

《 Disinfectable Surfaces 》

※Do not use this product on leather or artificial leather.
※This product may discolor the surfaces wiped.

[Plastic] Nylon, polypropylene, soft polyethylene, hard polyethylene, polyethylene, ABS, acrylic, polycarbonate, PET, PVC, epoxy [Metal] Stainless steel, aluminum
[Rubber] NBR, EPDM, butyl rubber, urethane, silicone [Other material] Glass

Cautions

●Do not flush wipes down the toilet. ●Close the lid after using. Use up the wipes quickly after opening. ●Avoid contact with eyes and mucous. ●Stop using this product if any syndrome occurs. ●Do not use this product on leather or artificial leather goods. ●Confirm the instruction manual before using the wipes on a machine. ●Wear gloves when using the product. ●Wash and dry the 250-count bucket once a month. ●Keep out of the reach of children. ●Avoid direct sunlight and wet environment. Store in a clean and cool place. ●The wipes may discolor after a long period of time from opening. The discoloration does not deteriorate the disinfectant efficacy.

Product Code	Specifications	No. of Wipes	Qty/carton	JAN Code
2600270	150×300mm	30 pcs	36 packs	4 958286 078397
2600271	150×300mm	60 pcs	18 packs	4 958286 078403
2600272	Bucket 150×300mm	250 pcs	4 buckets	4 958286 078410
2600273	Refill Pack 150×300mm	250 pcs	6 packs	4 958286 078427

Made in Japan



Design Your Voice
Hakuzo

Add: 2-4-9, Tokui-cho, Chuo-ku,
Osaka 540-0025, Japan

TEL: 06-6942-0451

FAX: 06-6945-0539

www.hakuzo.co.jp/eng/index.php