

Environment Disinfection Use Wet Gloth

HAKUZO





Environment Cloth Double Block

Non-woven fabric size $150_{mm} \times 300_{mm}$

30 sheets pack 60 sheets pack Now Available

Wet cloth for environment disinfection with US EPA certified Quaternary Ammonium Salt Impregnation.

- ●60pcs pack with one push cap.
- ■30pcs pack which can be hooked by S-hook.



Large Sheet Packet (60pcs)

Size: 150mm×300mm



Large sheet Container (250pcs)

Size: 150mm×300mm



Large Sheet Refill Packet (250pcs Size: 150mm×300mm





Can be hanged on rack or something similar.

Large Sheet Packet (30pcs)
Size: 150mm×300mm

Just one simple push!

EPA certified quaternary ammonium salt impregnated

Environment Disinfection Use Wet Cloth

HAKUZO





Environment Cloth Double Block









30 sheets pack

60 sheets pack

250 sheets Refill Pack

250 sheets pack

Size

150mm×300mm

Base Material
Content

Rayon Non-woven fabric

Quaternary Ammonium Salt, water

Product Features

- Solution approved by US EPA is being used, effective in getting rid of a wide range of bacteria and germs.
- Almost no effect on material wiped, environmental friendly.
- Solution is neutral (ph of around 7.7).



- Getting rid of bacteria and germs around table, doorknob, bed and etc.
- Cleaning and getting rid of bacteria and germs on contaminated areas







Warning and Precautions

- Please keep the cap properly closed after using. Once the seal is being torn off, please use up as soon as possible.
- •Please do not use directly on eyes, and mucous membranes.
- Please do not flush in the toilet bowl.
- This product is not intended to be used on cleaning body.
- •Please stop using immediately if symptoms persist.
- Please do not use on leather goods.
- Before using on machine, please read the cleaning manual for the machine.
- When refill, or around once a month, please confirm whether is the packet or containerand cap are in good condition.

StoragePrecautions

- ●Please keep it out of children's reach.
- ●Please keep out of hot and humid environment.

 Avoid direct sunlight and store in a clean environment.

| Product No. | Item | Size | No. of sheets | Qty/case | JAN Code |
|-------------|---------------------------|---------------|---------------|-----------------|-----------------|
| 2600170 | Large Sheet Packet | 150mm × 300mm | 30sheets | 36pack/case | 4 958286 064154 |
| 2600173 | Large Sheet Packet | 150mm × 300mm | 60sheets | 18pack/case | 4 958286 064161 |
| 2600174 | Large Sheet Packet | 150mm × 300mm | 250sheets | 4container/case | 4 958286 061795 |
| 2600175 | Large Sheet Refill Packet | 150mm × 300mm | 250sheets | 6container/case | 4 958286 061801 |



Add: 2-4-9, Tokuicho, Chuo-ku Osaka 540-0025, Japan

TEL: 06-6942-0451 FAX: 06-6945-0539

www.hakuzo.co.jp/eng/index.php