

Maximum Barrier Precaution Set

MBPキットI

Sterile

Sterile kit with all the required gears for CVC insertion.



Tray type



Sterile bag type


- CVC insertion comes with high infection risks. This kit can help to reduce infection risks.
- Custom made orders are welcome.
- Other than mask and cap, all the other gears are sterile.

Main components


- Cap
- Mask
- Isolation gown
- Drape with hole
- Band aid

MBP set operation manual


STEP.1 Wearing of mask and cap




STEP.2 Opening the sterile wrap



STEP.3 Putting on an isolation gown and gloves



STEP.4 Opening the drape



Hakuzo Medical Corporation

4-9, 2-chome, Tokuicho, Chuo-ku, Osaka
540-0025, Japan
Tel; 81-6-6809-4319 Fax; 81-6-6945-0539
E-mail; hmc-trade@hakuzo.co.jp
<http://www.hakuzo.co.jp/eng/>