

Ophthalmology absorbent sponge

HAKUZO EYE SPONGE


Eye sponge | Ophthalmic sponge

Sterile

Product Feature

Big absorbent part, high in absorbency

To absorb fluid during an eye operation, making the operation smoother and safer



Size and Ingredient

Absorbent part: About 8mm×16mm
Ingredient: Cellulose



Before absorbing



After absorbing



10 pcs × 40 packs

Usages

- To absorb blood, and other fluids.
- To press on conjunctiva sclera.

使用期限：製造日から5年間

Specifications

Product No.	Packaging	Size/weight	Qty/case	Size/weight	JAN Code
2600401	10pcs×40packs	W162×D102×H65mm	12Boxes	W274×D227×H327mm	4 958286 072791
		92g		1,400g	