



Design Your Voice
Hakuzo

New Release

Ophthalmology absorbent sponge


HAKUZO EYE STICK

EYE STICK | Ophthalmic sponge spear

Sterile

Product Feature ①

Big absorbent part, high in absorbency



To absorb fluid during an eye operation, making the operation smoother and safer

Size


Absorbent part: About 8mm x 16mm
Whole body: About 85mm

Product Feature ②

Stick is big, and hexagon in shape, making gripping easy and comfortable. No pincettes are required, can be used immediately after opening the packaging.



Tip shape: hexagon



Absorbent part: cellulose
Handle part: polypropylene


5 sticks x 20 packs



Usages

- To absorb blood, and other fluids.
- To press on conjunctiva sclera.

Specifications

Product No.	Packaging	Size/weight	Qty/case	Size/weight	JAN code
2600400	5sticks x 20packs	W155×D85×H105mm	12 boxes	W274×D227×H327mm	 4 958286 072203
		172g		2,335g	



Design Your Voice
Hakuzo

Add: 2-4-9, Tokuicho, Chuo-ku
Osaka 540-0025, Japan
TEL: 0081-6-6942-0451
FAX: 0081-6-6945-0539
www.hakuzo.co.jp/eng/index.php