

Medical Device




Cerebral Sheet XJ

Single Use

EO Sterile

Neurosurgical patties | X-ray detectable

Nonwoven pads for neurosurgery. Comes with x-ray detectable threads.




Product Features

- Long staple cotton is being used.
- Nearly lint-free, excellent in water retention and water absorption.
- EO sterile.
- Highly visible x-ray detectable thread is being used.

Purposes

- To protect brain tissue, absorb fluid, stop bleeding, etc. during neurosurgery.

Specifications



5 pads are stacked and fixed to the left and right of the paper holder.

Product Code	Sizes	Qty/box	Box Size	Weight/box	Qty/ Carton	Carton Size	Weight/carton	JAN Code
2640121	10x13mm 10pcs/pouch	20 pouches	W145mm	Abt 160g	12 boxes	W455mm	Abt 2.30kg	4958286081908
2640122	10x30mm 10pcs/pouch	20 pouches	x	Abt 165g	12 boxes	x	Abt 2.35kg	4958286081915
2640123	20x20mm 10pcs/pouch	20 pouches	D125mm	Abt 170g	12 boxes	D275mm	Abt 2.40kg	4958286081922
2640124	20x40mm 10pcs/pouch	20 pouches	x	Abt 185g	12 boxes	x	Abt 2.60kg	4958286081939
2640125	30x60mm 10pcs/pouch	20 pouches	H140mm	Abt 220g	12 boxes	H320mm	Abt 3.00kg	4958286081946



Design Your Voice

Hakuzo

Add: 2-4-9, Tokui-cho, Chuo-ku, Osaka 540-0025, Japan

TEL: 06-6942-0451

FAX: 06-6945-0539

www.hakuzo.co.jp/eng/index.php