


**Sterile**

# CEREBRAL SHEET

Neurosurgical Patties | Non-woven pads for neurological surgery



Five pieces are stacked each on left and right of the backing paper, wrapped and fixed.

## Product Features

- Made from 100% long fiber cotton that hardly fluffs. Excellent in water absorption and water retention.
- Easy to conform to surface tissue and easy to peel off.
- A wide range of sizes is available. Can be used for different purposes.


## Usage

- Protection of nerve tissue, absorption of fluid, hemostasis, etc. during surgery.

## Specifications

Product No.	Size and Packaging	Qty/box	Box Size	Box Weight	Qty/carton	Carton Size	Carton Weight	JAN Cartons
2640203	10mm × 10mm 10 pcs/sheet/pack	20 packs	W145mm × D125mm × H140mm	154 g	12 boxes	W455mm × D272mm × H317mm	2.35kg	4958286074528
2640204	10mm × 30mm 10 pcs/sheet/pack	20 packs		161 g	12 boxes		2.4kg	4958286074535
2640205	20mm × 20mm 10 pcs/sheet/pack	20 packs		165 g	12 boxes		2.45kg	4958286074542
2640206	20mm × 40mm 10 pcs/sheet/pack	20 packs		178 g	12 boxes		2.65kg	4958286074559
2640207	30mm × 60mm 10 pcs/sheet/pack	20 packs		214 g	12 boxes		3.05kg	4958286074566
* 2640201	5mm × 5mm 10 pcs/sheet/pack	20 packs		152 g	12 boxes		2.32kg	4958286074504
* 2640202	7mm × 7mm 10 pcs/sheet/pack	20 packs		153 g	12 boxes		2.34kg	4958286074511

\* Built to order Products



Design Your Voice

**Hakuzo**

Add: 2-4-9, Tokuicho, Chuo-ku,  
Osaka 540-0025, Japan

TEL: 06-6942-0451

FAX: 06-6945-0539

[www.hakuzo.co.jp/eng/index.php](http://www.hakuzo.co.jp/eng/index.php)