



LAPAROTOMY SPONGES

Cylindrical Sponge G

Cylindrical Sponge G can be used as laparotomy sponge, made from cellulose.

- Sterilized with AC
- Impregnated with saline solution
- Contains X-ray detectable thread



Size	Packing	Amount of saline
φ 80mm × 280mm	6pcs/box, 4 boxes/carton	200mL
φ 50mm × 280mm	12pcs/box, 4 boxes/carton	100mL



Rapatto G Square gauze sponge

Rapatto G is gauze sponge with X-ray detectable thread and string.

- Available in required sizes and number of ply of gauze
- Available in both sterile and non-sterile

Size	Packing (non-sterile)
30 × 30cm 4ply	50pcs/ bag
45 × 45cm 4ply	50pcs/ bag



Hakuzo Medical Corporation
4-9, 2 chome, Tokuicho, Chuo-ku, Osaka
540-0025, Japan
Tel; 81-6-6809-4319 Fax; 81-6-6945-0539
E-mail; trade@hakuzo.co.jp