

Surgical Gauze Swab

Tupfer Swab G

For Surgical Use

[Medical Device] Gauze sponge with X-ray detectable thread
Product Name: Hakuzo Sterile Gauze Ball G



[Product Features]

- Comes with X-ray detectable thread G.
 - Can be customized and provided as both non-sterile and sterile swabs.
 - The position of x-ray detectable thread can be changed.
- ※This product is built to order.



| Product Code | Type | Size | Packaging |
|--------------|-----------------------------|------------|------------------------|
| 2270204 | L•5 balls/pack/sterile | 8mm x 15mm | 5 balls x 10 packs/box |
| 2270205 | S•5 balls/pack/sterile | 5mm x 10mm | 5 balls x 10 packs/box |
| 2270014 | L•100 balls/bag/non-sterile | 8mm x 15mm | 100 ball/bag |
| 2270015 | S•100 balls/bag/non-sterile | 5mm x 10mm | 100 balls/bag |