

Compressed Pad for Blood Stanching

MADE IN JAPAN



SP Cotton sterile

- Can be used for blood stanching after dialysis
- Also available in non-sterile

Size	Pcs/bag	Bags/box	Boxes/carton
18mm × 25mm	8	50	10



Appakumen sterile

Large-sized compressed pad

- Can be used for blood stanching after angiography
- Also available in non-sterile

Size	Pcs/bag	Bags/box	Boxes/carton
35mm × 60mm	1	50	10



Tablet Cotton Ball non-sterile

- Can be used for blood stanching after dialysis
- Also available in sterile

Product	Size	Pcs/bag
No.6	φ 20	100
No.30	φ 28	100



Tablet Cotton Ball with silicone non-sterile

This is a trial product.

Silicone is attached to compressed pad, which helps to increase the effect of blood stanching.

- Can be used for blood stanching after dialysis
- Also available in sterile

Size	Pcs/bag
φ 20	100



Hakuzo Medical Corporation

4-9, 2 chome, Tokuicho, Chuo-ku, Osaka

540-0025, Japan

Tel; 81-6-6809-4319 Fax; 81-6-6945-0539

E-mail; trade@hakuzo.co.jp