

Adhesive Bandage LFF

LATEX
FREE

Elastic adhesive bandage



【Materials】

Polypropylene, polyurethane





【Product Features】

- Can be fixed without tape. The position and tightness is readjustable.
- Non-adherent to clothing and hair.
- Can be cut by hand. No scissor is needed.
- Latex and acrylic free.
- Skin-friendly breathable porous bandage.
- Hydrophobic.

【Precautions for Use】

- Dry the skin before using this product.
- Do not apply this product directly to wounds.
- Stop using when any negative effects like rash, skin redness and itching occur.
- Overtightening may cause obstruction for blood flow.

【Usage】 Fixation for gauzes and pad.

Product No.	Specifications	Packaging	Qty/case	JAN Code
3087825	2.5cm×4.5m	1roll/pack	36 boxes	 4 958286 712567
3087850	5.0cm×4.5m	12packs/box	16 boxes	 4 958286 712574