

Transparent Film Dressing **MEDICAL DEVICE**

ELEBAN FILM



S10

10 μ in thickness

Product	Size
No.50	2 inch \times 2.8 inch (50mm \times 70mm)
No.60	2.4 inch \times 3.6 inch (60mm \times 90mm)
No.80	3.2 inch \times 4.8 inch (80mm \times 120mm)
No.100	4 inch \times 6 inch (100mm \times 150mm)



Transparent Film Dressing **MEDICAL DEVICE**

ELEBAN FILM

S

Product	Size
No.50	2 inch \times 2.8 inch (50mm \times 70mm)
No.60	2.4 inch \times 3.6 inch (60mm \times 90mm)
No.80	3.2 inch \times 4.8 inch (80mm \times 120mm)
No.100	4 inch \times 6 inch (100mm \times 150mm)
No.120	4.8 inch \times 10 inch (120mm \times 250mm)



EIEBAN FILM S10 and ELEBAN FILM S are made of a thin film of vapor permeable polyurethane.

S10 and S are permeable to oxygen and moisture vapor, but impermeable to bacteria and fluid.

- Waterproof and bacteria-proof
- Excellent at elasticity and flexibility
- Packed individually, sterilized with EO
- Low-allergic adhesive
- Easy application and excellent at fixing

MADE IN JAPAN

Hakuzo Medical Corporation

4-9, 2 chome, Tokuicho, Chuo-ku, Osaka

540-0025, Japan

Tel; 81-6-6809-4319 Fax; 81-6-6945-0539

E-mail; trade@hakuzo.co.jp