



For baby's eyes protection during phototherapy treatment

Eleban Barrier

10 pcs pack

EO Sterile



To protect a baby's delicate eyes during phototherapy treatment



Product Features

[Materials]
97% cotton, acrylic adhesive, 3% polyurethane

- Aluminum foils are used to block the lights during treatment.
- Adhesive tape is used to fix the eye mask, it is easy to put on and take off.
- Cotton is the main material, with acrylic adhesive that is gentle to the skin is used.

Product code	Type	Specifications	Packaging	Qty/case	JAN code
1330201	Newborn	5cm×12.5cm(Pad 2.5cm×9cm)	10 pcs/box	10 boxes	 4 958286 510019
1330202	Premature infant	5cm×10.5cm(Pad 1.8cm×7.5cm)	10 pcs/box	10 boxes	 4 958286 510026



Hakuzo Medical Corporation
4-9, 2-chome, Tokuicho, Chuo-ku, Osaka
540-0025, Japan
Tel: 81-6-6809-4319 Fax: 81-6-6945-0539
E-mail: hmc-trade@hakuzo.co.jp
<http://www.hakuzo.co.jp/eng/>