

Hakuzo Phototherapy Eye Protector

Silicone Adhesive

Eye mask for
protecting infant's eyes during phototherapy.

EO
Sterile


【Materials】 Cellulose nonwoven fabric, aluminum foil, rayon nonwoven fabric, silicone adhesive, PET film.

【Product Features】

- Aluminum foil can protect eyes from ultraviolet light.
- Hypoallergenic silicone adhesive is being used.
- Adhesive eye mask can conform to infant's eye surrounding easily.
- Easy to adjust the position with the large size of aluminum foil.

【Sizes】

- 2 sizes are available.



Specifications

10 pcs

Product Code	Specifications	Qty/box	Qty/carton	JAN Code
1400070	Newborn 50mm × 125mm	10 pcs/box	15 boxes	4 958286 080024
1400071	Low birthweight newborn 50mm × 105mm	10 pcs/box	15 boxes	4 958286 080031

Made in Japan



Design Your Voice
Hakuzo

Add: 2-4-9, Tokui-cho, Chuo-ku,
Osaka 540-0025, Japan

TEL: 06-6942-0451

FAX: 06-6945-0539

www.hakuzo.co.jp/eng/index.php