

Excellent absorbency keeps
your skin dry and clean

SANGO PAD NATURAL

MATERNITY PAD

Superabsorbent
Polymer
Skin Friendly Fabric

Maternity pad with Super Absorbent Polymer.



① Skin Friendly

Spongy pad and smooth surface can protect your delicate skin.

② Superb absorbency

With the inclusion of superabsorbent Polymer, this product can absorb lochia immediately and keep your skin dry and clean.

③ Prevent leaking

Both sides come with waterproof sheets that help prevent leakage.

Specifications

Product No.	Product Name	Size	Qty/Weight/Absorbency	JAN Code
1400014	meetuu Sango Pad Natural	L(15cm×52cm)	1pc/pack×5	4 958286 709352
			about 42g/pc about 680g/pc	
1400015		M(12cm×32cm)	1pc/pack×10	4 958286 709369
			about 18g/pc about 440g/pc	
1400016		S(7.5cm×23cm)	1pc/pack×20	4 958286 709376
			about 11g/pc about 160g/pc	

※Normal saline solution is used to test the absorbency.



Design Your Voice

Hakuzo

Add: 2-4-9, Tokuicho, Chuo-ku,
Osaka 540-0025, Japan

TEL: 06-6942-0451

FAX: 06-6945-0539

www.hakuzo.co.jp/eng/index.php