

Can be used for checking, nursing, treatment purposes.

A convenient multipurpose sheet.

Now on sale



HAKUZO CURE CARE SHEET

High absorbency disposable sheet



Absorbing surface


Waterproof surface

【Product Features】


- Waterproof sheet is being used to prevent leakage.
- High absorbent polymer is being used to prevent flow back and to keep the surface dry.
- 3 sizes are available for different usages.

【Absorbency level】

※Saline was used for this experiment



【Product Structure】



【Main Usages】

Cure (Trmenteat)

- For use during surgery
- For use during endoscopic examination
- For use during wound cleaning
- For use during delivery
- For use on newborn baby immediately after delivery

Care(Nursing)

- For use during the cleaning of genital area and buttocks
- For use of drying hair after shower.
- For use when changing diapers.
- A solution for incontinence and bedwetting.
- A solution for vomit

Product Code	Type	Specificaition	Packaging	Qty per case	JAN Code
3147001	S size	32cm × 44cm	200pcs/bag	4bags/case	4 958286 706870
3147002	M size	44cm × 59cm	100pcs/bag	4bags/case	4 958286 706887
3147003	L size	60cm × 90cm	40pcs/bag	5bags/case	4 958286 706894



Hakuzo Medical Corporation

4-9, 2-chome, Tokuicho, Chuo-ku, Osaka
540-0025, Japan

Tel: 81-6-6809-4319 Fax: 81-6-6945-0539

E-mail: hmc-trade@hakuzo.co.jp

http://www.hakuzo.co.jp/eng/