



A convenient multipurpose sheet.



HAKUZO CURE CARE SHEET

High absorbency disposable sheet



【Absorbency level】

※Saline was used for this experiment

S Size (30cm x 45cm) About 250mL	
M Size (45 cm x 60 cm) About 570mL	
L Size (60 cm x 90 cm) About 1600mL	

【Product Structure】

Absorbing surface	
	Hydrophilic spunbond nonwoven
	Absorbing paper
	Super Absorbent Polymer (SAP)
	Fluff pulp
	Absorbing paper
	Polyethylene film
Waterproof surface	

【Product Features】

- Waterproof sheet is being used to prevent leakage.
- High absorbent polymer is being used to prevent flow back and to keep the surface dry.
- 3 sizes are available for different usages.

Product Code	Type	Specificaion	Packaging	Qty per case	JAN Code
3147015	S size	30cm x 45cm	200pcs/bag	4bags/case	 4 958286 721057
3147016	M size	45cm x 60cm	100pcs/bag	4bags/case	 4 958286 721064
3147017	L size	60cm x 90cm	40pcs/bag	4bags/case	 4 958286 721071



Hakuzo Medical Corporation

4-9, 2-chome, Tokuicho, Chuo-ku, Osaka
540-0025, Japan
Tel: 81-6-6809-4319 Fax: 81-6-6945-0539
E-mail: hmc-trade@hakuzo.co.jp
<http://www.hakuzo.co.jp/eng/>