


Thin absorbent disposable sheets for multiple purposes



HAKUZO Absorbent Disposable Sheets

Medical Disposable Sheets



Product code	Size	Packaging	Qty/carton	JAN Code
3096000	45cm×60cm	50pcs/bag	4 bags	4 958286 712970
3096001	60cm×90cm	20pcs/bag	4 bags	4 958286 712987
3096002	90cm×90cm	10pcs/bag	6 bags	4 958286 712994
3096003	90cm×120cm	10pcs/bag	4 bags	4 958286 713007
3096004	90cm×200cm	10pcs/bag	3 bags	4 958286 713014

Size	Absorbent volume (saline)
about 45cm×60cm	about 1,200mL
about 60cm×90cm	about 2,300mL
about 90cm×90cm	about 3,500mL
about 90cm×120cm	about 4,500mL
about 90cm×200cm	about 8,200mL

- 【Product Features】**
- Comes with waterproof film to prevent leakage.
 - Highly absorbent polymer keeps the nonwoven surface dry.
 - Comes with 5 sizes, can be used for different purposes.



Add: 2-4-9, Tokui-Cho, Chuo-ku, Osaka 540-0025
 TEL: +816-6809-4319
 FAX: +816-6942-0496
 www.hakuzo.co.jp