

Light

Suction

Care


3in1

Refreshing

# Suction Oral Care Set

With Pen Light 

Multi-Functional Set




Slightly bent tip can fit gum firmly!

Can use control port to adjust the suction with 1 finger.

Can connect to tubes with diameter from 5~7mm.  
※Suction tube is not included.

Can illuminate the oral cavity!



IT IS SUITABLE FOR 

**WHO WANTS TO SAVE TIME**  
Multi-functional oral care set can improve the efficiency.

**WHO WANTS BETTER VISIBILITY AND ACCESS TO THE ORAL CAVITY**  
By illuminating the oral cavity, caregivers can have better visibility and access to the oral cavity.

**WHO WANTS TO REDUCE THE RISK OF ASPIRATION WHEN GIVING ORAL CARE**  
By connecting the swab/toothbrush to a suction equipment, the contaminants can be removed efficiently and the risk of aspiration can be largely reduced.


Multi-functional!

Comes with a pen light!

Single hand operable!

point  
1

## Ultra Bright Pen Light



**Can illuminate the whole oral cavity!**

With the use of ultra bright pen light, the oral cavity can be fully illuminated to avoid overlooking any debris and bleeding spots.



point  
2

## Soft Bristles




**Can be used to clean mucosa and gums.**

**0.076mm** ultra fine bristles are being used. Can be used to clean mucosa and gums. Small tip can access hard-to-reach areas in the mouth easily.


point  
3

## Small Sponge Tip



**Can access hard-to-reach areas in the mouth.**

**Small and round shaped sponge tip** is suitable for patients with trismus. **Jagged and firm sponge** is efficient in removing debris.

Product Name	Product Code	Qty/Box	Qty/Carton	JAN Code
Refreshing Suction Oral Care Set with Pen Light	1850830	6 sets	10 boxes	 4 958286 720012



*Design Your Voice*

# Hakuzo

Add: 2-4-9, Tokui-cho, Chuo-ku,  
Osaka 540-0025, Japan

TEL: 06-6942-0451

FAX: 06-6945-0539

[www.hakuzo.co.jp/eng/index.php](http://www.hakuzo.co.jp/eng/index.php)