



Cotton tipped swab for oral care

# Soukai Menbou

Cotton tipped swab available in plural care




## Features

**Sponge**  
Jagged-cut shape, which is easy to care for back of tongue or coat of tongue

**Shaft**  
Uses thick and strong paper

**Cotton Tip**  
Drum type, which is easy to care for a wide range of sites as in buccal mucosa

## How to use



Packing
5 pcs/bag, 20 bags/box, 10 boxes/carton

**MADE IN JAPAN**



**Hakuzo Medical Corporation**  
 4-9, 2 chome, Tokuicho, Chuo-ku, Osaka  
 540-0025, Japan  
 Tel; 81-6-6809-4319 Fax; 81-6-6945-0539  
 E-mail; trade@hakuzo.co.jp

