

Assembled in Japan

Skin Friendly One Touch Binder

To protect, warm
and fix abdomen after surgery.

Bleached Cotton Cloth


To fix stoma

Abdominal
Binder
for
Postoperative
Care


Easy and comfortable

It fits body well and
the fabric is skin friendly.

Superior Stretchability




Easy to wear,
excellent in
breathability.




Bleached Cotton Cloth

Skin Friendly One Touch Binder



【The front】



Product Features

- Breathable fabric barely causes skin rash.
- Fits body well and can be fixed firmly.
- With suitable stretchability, the binder can offer mild compression on the wound.

Instructions

- Confirm the front and back sides before putting on. (Refer to the image above)
- The bleached cotton cloth should be facing outward.



Purposes

- To protect, warm and fix abdomen after surgery.
- Gauze fixation
- Stoma fixation
- For other abdominal protection.

Materials

Main material: Cotton
Stretchable material: Rayon, polyurethane



Product Code	Size	Suitable Waist Sizw	Packaging	Qty/carton	JAN Code
1730130	Regular (31×140cm)	75cm~100cm	1 pc/pack	30 packs/carton	 4 958286 715230
1730131	Small (31×120cm)	55cm~80cm	1 pc/pack	30 packs/carton	 4 958286 715247



Design Your Voice

HAKUZO

Add: 2-4-9, Tokui-Cho, Chuo-ku,
Osaka 540-0025

TEL: +816-6809-4319

FAX: +816-6942-0496

www.hakuzo.co.jp